

















## *Arabidopsis thaliana* Culture in Hydroponics

### Principle

Culture of *Arabidopsis thaliana* on a hydroponics system enabling experiments where roots can be harvested relatively easily and without excessive damage.

### Hoagland Solution and variants

#### Equipment and reagents

Machine/Product	Reference (Company, Type, ...)
Volumetric flasks	250 ml, 1 L, 2L
Hotplate/Stirrer	 Jenway Model 1100
Potassium Nitrate KNO <sub>3</sub>	 Merck 60410
Calcium Nitrate Ca(NO <sub>3</sub> ) <sub>2</sub> •4H <sub>2</sub> O	 Vel 6560
Ammonium Phosphate NH <sub>4</sub> H <sub>2</sub> PO <sub>4</sub>	 Merck 1.01126
Magnesium Sulfate MgSO <sub>4</sub> •7H <sub>2</sub> O	 Merck 105886
Ferrous Sulfate FeSO <sub>4</sub> •7H <sub>2</sub> O	 Merck 1.03965
Na <sub>2</sub> -EDTA (disodium salt dehydrate) 99 %; 372.2 g.mol <sup>-1</sup>	 Sigma ED2SS
Boric Acid H <sub>3</sub> BO <sub>3</sub>	 Baker 0501
Manganese chloride MnCl <sub>2</sub> •4H <sub>2</sub> O	  Merck 1.05927
Copper sulfate CuSO <sub>4</sub> •5H <sub>2</sub> O	  Baker 0104
Molybdcic Acid H <sub>2</sub> MoO <sub>4</sub>	 Aldrich 232084
Zinc sulfate ZnSO <sub>4</sub> •7H <sub>2</sub> O	   UCB1884

#### Hoagland Stock Solutions

Weigh off the different reagents for each stock as given in table below:

Macro elements			1 L /10 L	H 1:1
	g.mol <sup>-1</sup>	g / 2 L		mM
KNO <sub>3</sub>	101,11	20.4		100
Ca(NO <sub>3</sub> ) <sub>2</sub> •4H <sub>2</sub> O	236,15	14.16		30
MgSO <sub>4</sub> •7H <sub>2</sub> O	246,48	9.8		20

<b>Phosphorus</b>			1L / 10L	
NH <sub>4</sub> H <sub>2</sub> PO <sub>4</sub>	115,03	4.6		20
<b>Iron Solution</b>		g / 250 ml	6ml / 10L	
FeSO <sub>4</sub> •7H <sub>2</sub> O	278,02	1,9		27
Na <sub>2</sub> -EDTA•2H <sub>2</sub> O	372,2	1,25		13
<b>Micro elements</b>		g / 1 L	10 ml /10 L	
H <sub>3</sub> BO <sub>3</sub>	61,83	2,86		46
MnCl <sub>2</sub> •4H <sub>2</sub> O	197,91	1,81		9.1
CuSO <sub>4</sub> •5H <sub>2</sub> O	249,68	0,08		0.32
H <sub>2</sub> MoO <sub>4</sub>	161,97	0,09		0.55
ZnSO <sub>4</sub> •7H <sub>2</sub> O	287,54	0,22		0.76

**Note:** Heating can be necessary to dissolve molybdate in the micronutrients solution. Work on a laboratory hotplate with constant stirring.

**Note:** Stock solutions are autoclaved after being prepared to prevent growth of algae and microfauna.

**Note:** Macro-elements stock solution tends to precipitate when staying at 4°C for longer periods. The precipitate can be dissolved by bringing to room temperature and putting in a ultrasonic bath for a few minutes.

#### *Hoagland Work Solutions*

##### **Hoagland 1:10**

Macro elements without phosphorus..... 100 ml for 10 L  
 Phosphorus solution..... 100 ml for 10 L  
 Iron Solution..... 0.6 ml for 10 L  
 Micro elements..... 1 ml for 10 L


##### **“HP” Hoagland**

Macro elements without phosphorus..... 100 ml for 10 L  
 Phosphorus solution..... 50 ml for 10 L  
 Iron Solution..... 0.6 ml for 10 L  
 Micro elements..... 1 ml for 10 L

##### **“LP” Hoagland**

Macro elements without phosphorus..... 100 ml for 10 L  
 Phosphorus solution..... 12.5 ml for 10 L  
 Iron Solution..... 0.6 ml for 10 L  
 Micro elements..... 1 ml for 10 L

#### **Hydroponic device set-up**

<b>Machine/Product</b>	<b>Reference (Company, Type, ...)</b>
Eppendorf Microtubes	1.5 ml
PE containers	Fisher Bioblock
Fitting, Perforated Lids	
Agar No.2 "Bacteriological"	 Lab M Ltd. MC006
Paper tape	


Autoclave	Tuttnauer
Hoagland Solutions	See 0 <i>Hoagland Work Solutions</i> above
Peristaltic Pump	
Climate Chamber	Snijders

### *Pre-treatment of hydroponics containers*

Plastic containers are made opaque. First they are coated with primer spray paint (Motip Primer Grey Eriks+Baudouin Catalog#04054) and the with black spray paint (Motip 4006 RAL9005 mat black; Eriks+Baudouin Catalog#04006).

When used for the first time the containers are prewashed and rinsed with 10 % HCl, followed by deionized water, to remove every wall bound metal to prevent metal contamination in the nutrient solution.

Fitting lids are made in opaque gray plastic (PE) with holes of approx. 11 mm diameter.




*Figure 1: Example of Hydroponics container and lid*

*Table 1: Types of containers mostly used for experiments with their nominal volume, their useful volume(max volume with lid + microtubes) and number of holes that could be fitted in the corresponding lids.*

Container type	Useful volume	# holes in lid
3 L	2900 ml	9x9 = 81
1.5 L	1400 ml	9x4 = 36
1 L	960 ml	4x4 = 16

### *Agar Seed Bed*



*Figure 2: Making and filling the agar seed bed devices*

Cut off the lid and the bottom half of standard 1.5 ml Eppendorf tubes

Autoclave the tubes at 121 °C

Place the tubes upside down on a strip of paper adhesive tape

Prepare the agar solution from “LP Hoagland” with 6 g Agar No2 per liter (0.6 %)
   
Autoclave the solution and leave it to cool to about 40-45 °C






Fill the tubes with agar solution and make sure they are filled over the brim “with a bubble” as pictured in Figure 2 above. Use a pipette or a

Let cool down to room temperature ( $\pm 15$  min) and transfer to the cut-out lids of the hydroponics containers.

Transfer the lids to the hydroponics containers filled with the appropriate pre-culture solution and place in a high-humidity environment (65 % RH).

### Arabidopsis thaliana culture

Machine/Product	Reference (Company, Type, ...)
Disposable Petri Dish	90 mM diameter
Filter Paper	S&S blue ribbon 5893
Fine-tipped tweezers	
Small diameter tubing	
Peristaltic Pump	
Tubing for Peristaltic Pump	
Climate Chamber	Snijders
Commercial Bleach	 LODA Javel <5% Cl 10°
Nitric Acid, concentrated HNO <sub>3</sub> (65 %; 1,41 g.ml <sup>-1</sup> ; 63.01 →14.55 mol.l <sup>-1</sup> )	  Merck 20429.320

### Vernalization (cold treatment)

Cold shock is a germination improvement technique where an incubation period of at 4 °C is applied to the seeds to synchronize germination.

For *Arabidopsis thaliana* the seeds are placed on wet filter paper in a closed Petri dish for 3 days at 4 °C in the dark before sowing.

### Sowing, Germination and Climatic Growth conditions

The seeds are transferred in pairs to the agar with fine tweezers seed beds placed in nutrient solution. The germination conditions are 14/10 day/night cycle at 165  $\mu\text{mol.m}^{-2}$ , with respectively 22 °C / 18 °C and a constant relative humidity of 65 %. (see Table 2 below)

It takes about 3 days for the seeds to germinate and after about 1 week the agar seed beds can be “thinned out”, which means only one seedling is left per tube and the excess is removed using fine tweezers.

Table 2: Climate chamber program for *Arabidopsis thaliana*

Process step				ChambE	HmdtyE	illumE
ChambS	HmdtyS	illumS	Time			
<b>1. All but 24</b>				22 °C	65 %	65 %
22 °C	65 %	65 %	05:30			15:00 → 20:30
<b>2. All but 24</b>				18 °C	65 %	10 %
22 °C	65 %	65 %	00:30			20:30 → 21:00
<b>3. Dark</b>				18 °C	65 %	10 %
18 °C	65 %	10 %	09:30			21:00 → 06:30

<b>4. All but 24</b>				22 °C	65 %	65 %	06:30 → 07:00
18 °C	65 %	10 %	00:30				
<b>5. All but 24</b>				22 °C	65 %	65 %	07:00 → 15:00
22 °C	65 %	65 %	08:00				
<b>6. Repeat</b>				0.0	0.0	0.0	
1.0	0.0	99.0	00:00				

### *Maintenance of the Hydroponics Culture*

Once the Arabidopsis plants are a week old, the nutrient solution has to be replaced twice a week. For experimental pre-culture, "Hoagland HP" is usually used.

From the time the plants are 1 week old, the nutrient solution also has to be aerated. This is achieved by placing a piece of small-diameter diameter about 5 cm in the nutrient solution and bubbling air by means of a peristaltic pump placed next to the growth cabinet.

### *Dismantling and cleaning the hydroponics setup.*

#### **Tubes:**

- Remove the agar from the Eppendorf tubes. Place the tubes in a 2 L PE bottle and fill with a detergent solution and a dash of household bleach. Shake overnight.
- Rinse the tubes with hot tap water, then with demineralized water. Put the clean tubes in a glass beaker and steam-sterilize them in the autoclave.

#### **Containers and Lids:**

- Remove all labels and writings from the containers and lids. (Do not use acetone as it corrodes the plastic from the lids). Wash them thoroughly in hot water and detergent. Make sure to scrub the corners and the bottom since precipitates accumulate there most easily.
- [optional] If the experiment carried out contained heavy metals or radionuclides, the containers have to be rinsed with acid solution ( $\text{HNO}_3$  1M = 70 ml /L). Use a shallow platter for this.
- Rinse once with tap water and rinse once with dilute household bleach (<5%). Use a shallow platter for this.
- Thoroughly rinse with demineralized water to remove any trace of bleach.